



## Model EMPX-HF ELECTROMAGNET-BASED HORIZONTAL FIELD CRYOGENIC PROBE STATION

### Introduction

The Model EMPX-HF is a versatile cryogenic electromagnet-based micro-manipulated probe station used for non-destructive testing of devices on full and partial wafers up to 25 mm (1 in) in diameter. The EMPX-HF is a platform for measurement of magneto-transport, electrical, electro-optical, parametric, high Z, DC, RF, and microwave properties of materials and test devices. Nanoscale electronics, quantum wires and dots, semiconductors, superconductors, and spintronic devices are typical materials measured in an EMPX-HF. A wide selection of probes, cables, sample holders, and options makes it possible to configure the EMPX-HF to meet your specific measurement applications.

Mounted on a vibration isolation table, the EMPX-HF is equipped with a 5.5 kOe (0.55 T) horizontal (in-plane) field electromagnet. The EMPX-HF operates over a temperature range of 4.5 K to 400 K. With options, the base temperature can be extended down to 3.2 K. The probe station provides efficient temperature operation and control with a continuous refrigeration system using either liquid helium or liquid nitrogen. Field dependent measurements at ambient temperature are possible without the use of cryogens. Liquid helium is only required for cooling the device under test (DUT) to below 80 K. Vapor-cooled shielding optimizes efficiency and intercepts blackbody radiation before it reaches the sample. A heater on the sample stage along with independent control of radiation shield heaters provides the probe station with fast thermal response.

The EMPX-HF is user configured with up to four ultra-stable micro-manipulated stages, each providing precise 3-axis control of the probe position to accurately land the probe tip on device features. The 360° sample stage rotation option allows you to measure angular-dependent and anisotropic magneto-transport properties of the DUT. Proprietary probe tips in a variety of sizes and materials minimize thermal mass and optimize electrical contacts to the DUT. Probe tips are thermally linked to the sample stage to minimize heat transfer to the DUT.

For increased versatility, EMPX-HF options include a 3.2 K base temperature stage, 360° sample stage rotation, a LN<sub>2</sub> Dewar kit, higher magnification microscope, vacuum turbo pumping system, pump-line vibration isolator, recirculating chillers, and fiber optic probe arm modification.

### Features

- 5.5 kOe (0.55 T) horizontal (in-plane) field electromagnet
- 360° sample stage rotation option
- High stability operation from 3.2 K to 400 K
- Measurements from DC to 67 GHz
- Accommodates up to 25 mm (1 in) diameter wafers
- Configurable with up to four thermally anchored micro-manipulated probe arms
- Probe arms with 3-axis adjustments and  $\pm 5^\circ$  planarization
- Cables, shields, and guards minimize electrical noise and thermal radiation losses
- Options and accessories for customization to specific research needs

#### Lake Shore Cryotronics, Inc.

575 McCorkle Boulevard ■ Westerville, Ohio 43082

Phone 614-891-2244 ■ Fax 614-818-1600

[info@lakeshore.com](mailto:info@lakeshore.com) ■ [www.lakeshore.com](http://www.lakeshore.com)

*All specifications subject to change. Rev 7/7/2010*

## Applications

- Angular-dependent and anisotropic magneto-transport measurements
- Electrical and electro-optical measurements over a wide temperature range
- RF and microwave
- Parametric testing
- Shielded/guarded/low noise characterization
- High Z
- Non-destructive, full wafer testing

## Materials

- Nanoscale electronics (carbon nanotube transistors, single electron transistors, molecular electronics, nanowires, etc.)
- Quantum wires and dots, quantum tunneling
- Single electron tunneling (Coulomb blockade)
- Superconductors
- Spintronic devices
- Basic semiconductor devices including organics, LEDs, and dilute magnetic semiconductors

## Specifications

### Magnetic Field

Maximum field	5.5 kOe (0.55 T)*
Probe movement due to magnet field ramping to 5.5 kOe (0.55 T)	<5 μm
<b>Field uniformity (at 70 A)</b>	
10 mm diameter wafer	0.6%
12.7 mm diameter wafer	1.1%
20 mm diameter wafer	2%
25.4 mm diameter wafer	2.6%
<b>Field control stability</b>	<50 mG RMS**

\*Hall probe is calibrated to read field at the center of the magnet gap

\*\*With settle band of 100 mG

### Temperature

Sample temperature range—standard	4.5 K to 400 K*	
Sample temperature range—optional	3.2 K to 400 K**	
<b>Temperature control (heaters)</b>		
Sample stage	100 W	
Radiation shields	Two, 100 W each	
Probe arm	Measurement only	
<b>Temperature control stability</b>		
With LHe	5.5 K to 50 K	<5 mK RMS
	50 K and higher	<50 mK RMS
With LN <sub>2</sub>	80 K and higher	<100 mK RMS

\*Limited from 8 K to 400 K when configured with the PS-360-EMPX 360°

sample stage rotation option


\*\*Selectable equipment

### Probe Arm Adjustments

<b>Travel</b>	
X axis	51 mm (2 in)
Y axis	25 mm (1 in)
Z axis	18 mm (0.7 in)
<b>Translation scale</b>	
X axis	20 μm
Y and Z axes	10 μm
Planarization*	±5°
<b>Sample stage (sample holder) in-plane rotation**</b>	
	360°

\*Included with microwave probes

\*\*Optional — recommended when taking measurements using microwave probes



**NOTE:** because of the location of the magnet poles in the EMPX-HF, the probes are aligned at 30° angles, which limits probe tip proximity (probe-dependent). Contact Lake Shore for more information.

## Frequency Range

<b>ZN50 DC/RF probe frequency range</b>	
<b>Tungsten with cryogenic coaxial cable</b>	0 to 50 MHz*
<b>Tungsten with semirigid coaxial cable</b>	0 to 1 GHz*†
<b>Paliney 7 with cryogenic coaxial cable</b>	0 to 50 MHz*
<b>Paliney 7 with semirigid coaxial cable</b>	0 to 1 GHz*†
<b>BeCu with cryogenic coaxial cable</b>	0 to 50 MHz*
<b>BeCu with semirigid coaxial cable</b>	0 to 1 GHz*†
<b>GSG microwave probe frequency range</b>	
<b>Low frequency with K connector</b>	0 to 40 GHz*
<b>Mid frequency with 2.4 mm connector</b>	0 to 50 GHz*
<b>High Frequency with 1.85 mm connector</b>	0 to 67 GHz*

\*Selectable equipment

† S21 > -10 dB up to 1 GHz, except for a (-40 dB) spike between 400 MHz and 800 MHz depending on probe model and placement; S11 < -3 dB up to 1 GHz

## Optical

<b>Optical viewport—located on top lids</b>	Ø54 mm (2.13 in) outer window and Ø50 mm (2 in) inner window
<b>Outer, clear fused quartz</b>	99% IR transmittance
<b>Inner</b>	IR absorbing with narrow band visible light transmittance
<b>Optical resolution—microscope</b>	
<b>7:1 zoom</b>	8 µm
<b>16:1 zoom</b>	8 µm*

\*Selectable equipment

## Sample Holder (Chuck)

<b>Maximum sample size—overall</b>	Up to Ø25 mm (1 in)
<b>SH-1.00-G-EM, grounded sample holder</b>	Up to Ø25 mm (1 in) and 400 K
<b>SH-1.00-I-EM, isolated sample holder</b>	Up to Ø25 mm (1 in) and 400 K*
<b>SH-1.00-C-EM, coaxial sample holder**</b>	Up to Ø25 mm (1 in) and 400 K*
<b>SH-1.00-T-EM, triaxial sample holder**</b>	Up to Ø25 mm (1 in) and 400 K*

\*Selectable equipment

\*\*Not available in combination with the 360° sample stage rotation feature

## Standard Equipment

<b>Electromagnet</b>	5.5 kOe (0.55 T), horizontal field
<b>Water flow rate</b>	7.6 L/min (2 gal/min)
<b>Electromagnet power supply</b>	Lake Shore Model 642
<b>Output type</b>	Bipolar, linear, 4-quadrant, DC current source
<b>Current</b>	±70 A
<b>Voltage</b>	±35 V
<b>Water flow rate</b>	5.7 L/min (1.5 gal/min)
<b>Gaussmeter and Hall probe</b>	Lake Shore Model 475 with probe (control and monitor field)
<b>Flow cryostat</b>	4.5 K to 400 K
<b>Sample stage temperature sensor</b>	Lake Shore Model CX-1050-SD-HT-1.4M calibrated Cernox™ RTD
<b>Sample stage heater</b>	100 W
<b>Cooled radiation shield and cooled IR-absorbing window above the sample</b>	
<b>Radiation shield temperature sensors</b>	Two Lake Shore Model DT-670C-CU silicon diodes
<b>Radiation shield heaters</b>	Two 100 W
<b>Removable top lid with viewport</b>	Ø50 mm (2 in) window
<b>Temperature control</b>	One Lake Shore Model 336 temperature controller and one Model 142 200 W (two channels, 100 W each) power supply (independent regulation of sample stage and radiation shield, and probe arm temperature monitoring)
<b>Electroless nickel-plated aluminum vacuum chamber</b>	
<b>Outside dimensions</b>	340 mm × 203 mm (13.4 in × 8 in)
<b>Removable top lid with clear fused quartz viewport</b>	Ø54 mm (2.13 in) window
<b>Probe ports</b>	4 surround the sample thermal radiation shield
<b>Pump port</b>	NW 40 (pump sold separately)
<b>Gas purge and 0.5 psi safety pop-off port</b>	NW 25
<b>Machined aluminum base plate</b>	787 mm × 419 mm (31 in × 16.5 in)
<b>PS-TMC vibration isolation system</b>	
	Includes pneumatically driven gimbal piston isolator, actuators, and supports; self leveling with a resonant frequency below 2 Hz; requires 40 psi nitrogen or air
<b>Isolator natural frequency</b>	Vertical 0.8 Hz, horizontal 1.0 Hz
<b>Isolator efficiency at 5 Hz</b>	Vertical 80 to 97%, horizontal 60 to 90%
<b>Isolator efficiency at 10 Hz</b>	Vertical 90 to 99%, horizontal 70 to 95%
<b>Temperature sensor installed and wired to a 6-pin feedthrough (included on one probe arm)</b>	
<b>Grounded sample holder</b>	SH-1.00-G-EM, accommodates up to a Ø25 mm (1 in) sample with a Ø25 mm (1 in) probe area
<b>Optics</b>	
<b>Zoom 70 microscope</b>	7:1 zoom with 8 µm resolution
<b>Color CCD camera</b>	S-video or composite output format
<b>Swing arm</b>	Optics can be manipulated to view any part of the sample or wafer, and can be retracted and swung away to allow access to the top of the vacuum chamber for sample exchange
<b>Video monitor</b>	High resolution, 17-inch
<b>Sample illumination</b>	Coaxial via fiber optic or ring light from an adjustable light source and power supply (must specify sample illumination at time of order) <b>NOTE: Coaxial illumination is recommended for highly reflective materials</b>
<b>High efficiency helium transfer line with foot valve for precise flow regulation</b>	
<b>Instrument console</b>	
<b>Basic tools and spares kit for standard operation</b>	

## Required User Configurable Equipment — Micro-manipulated Stages, Probes, Probe Tips, and Cables

### UP TO FOUR XYZ PRECISION MICRO-MANIPULATED STAGES

Part Number	Description
MMS-10	Micro-manipulated stage with thermal radiation shields, stainless steel welded bellows, and feedthrough ports—includes probe arm and base; probes, probe tips, and cables sold separately

### ZN50 DC/RF PROBES

- Ideal for: DC biasing, low/high frequency measurements, low noise shielded, and low-leakage guarded measurement
- ZN50 probe base incorporates a pair of copper braids that anchor to the sample stage to dynamically cool/heat the probe to the sample temperature
- SMA connector mounted directly to an alumina ceramic blade with a 50 Ω stripline routed to the probe contact

We understand that today's researcher requires flexibility. Our wide selection of probes, cables, sample holders, and options make it possible to configure a probe station to meet your specific measurement applications.

Part number (probe body)	Description
ZN50	50 Ω stripline probe body mount (each probe body mount requires a ceramic blade—selectable below)

Part number (probe body and ceramic blade)	Tip material	Maximum frequency (GHz)	Maximum probe temperature*	Maximum sample temperature**	Tip radius (μm)
ZN50R-03-W	Tungsten	1 Maximum frequency 50 MHz with ZN50C-G or ZN50C-T cable; maximum frequency 1 GHz with MWC-10-00K-NM cable	350 K	400 K	3
ZN50R-10-W					10
ZN50R-25-W					25
ZN50R-03-P7	Paliney 7				3
ZN50R-10-P7					10
ZN50R-25-P7					25
ZN50R-03-BECU	BeCu				3
ZN50R-10-BECU					10
ZN50R-25-BECU					25
ZN50R-100-BECU					100
ZN50R-200-BECU		200			

\* As measured by the probe arm temperature sensor

\*\* Selectable equipment

### ZN50 DC/RF CABLES

Part number	Cable type	Connector type	Feedthrough type	Measurement configuration	Maximum frequency	Maximum cable temperature*	Maximum sample temperature**
ZN50C-G	Ultra-miniature cryogenic coaxial	SMA	BNC	Shielded	50 MHz	350 K	400 K
ZN50C-T	Ultra-miniature cryogenic coaxial	SMA	3-lug triaxial	Low leakage	50 MHz		
MWC-10-00K-NM	Non-magnetic semirigid microwave coaxial	K (SMA compatible)	Loss-less compression seal	High frequency	1 GHz <sup>†</sup>		

\* As measured by the probe arm temperature sensor

\*\* Selectable equipment

<sup>†</sup> S21 > -10 dB up to 1 GHz, except for a (-40 dB) spike between 400 MHz and 800 MHz depending on probe model and placement; S11 < -3 dB up to 1 GHz

## GSG MICROWAVE PROBES\*

- Coplanar waveguide probe with ground-signal-ground (GSG) contact geometry
- User-specified pitch (spacing)
- Optimized low thermal conductivity coaxial leading to low thermal conductivity tips
- Cooled to near-sample temperature
- Separate planarization module with  $\pm 5^\circ$  rotation mechanism is provided

Part number	Connector type	Maximum frequency (GHz)	Maximum probe temperature* <sup>†</sup>	Maximum sample temperature**	Pitch (μm)
GSG-050-40A-E-NM	K	40	350 K	400 K	50
GSG-100-40A-E-NM					100
GSG-150-40A-E-NM					150
GSG-200-40A-E-NM					200
GSG-250-40A-E-NM					250
GSG-050-50A-E-NM	2.4 mm	50			50
GSG-100-50A-E-NM					100
GSG-150-50A-E-NM					150
GSG-200-50A-E-NM					200
GSG-250-50A-E-NM					250
GSG-050-67A-E-NM	1.85 mm	67			50
GSG-100-67A-E-NM					100
GSG-150-67A-E-NM					150
GSG-200-67A-E-NM					200
GSG-250-67A-E-NM					250

\* As measured by the probe arm temperature sensor

\*\* Selectable equipment

<sup>†</sup> PS-360-EMPX 360° sample stage rotation option; recommended when taking measurements using microwave probes

## GSG MICROWAVE CABLES

- Loss-less compression seal
- Semirigid with Teflon® dielectric

Part number	Cable type	Feedthrough type	Maximum cable temperature*	Maximum sample temperature**	Connector type	Maximum frequency
MWC-10-00K-NM	Non-magnetic semirigid microwave coaxial	Loss-less compression seal	350 K	400 K	K (SMA compatible)	40 GHz
MWC-10-240-NM					2.4 mm	50 GHz
MWC-10-185-NM					1.85 mm	67 GHz

\* As measured by the probe arm temperature sensor

\*\* Selectable equipment


## Sample Holders (Chucks)

Typical sample holder configuration characterized by:

- Leakage resistance between
  - Top surface and guard
  - Guard and ground
- Capacitance between
  - Top surface and guard
  - Guard and ground

### Types of sample holders

- Grounded sample holder—sample mount surface at system ground
- Isolated sample holder—backside contact not needed; sample mount surface is electrically non-conductive and isolated from ground
- Coaxial sample holder—backside contact can be made; sample mount surface is isolated from ground
- Triaxial sample holder—guarded backside contact can be made; sample mount surface has guarded isolation from ground



Part number	Measurement configuration	Separate feedthrough required	Maximum sample (diameter)	Maximum temperature
SH-1.00-G-EM	Grounded	No	Ø25 mm (1 in)	400 K
SH-1.00-I-EM	Isolated			
SH-1.00-C-EM ***	Coaxial	Yes*		
SH-1.00-T-EM ***	Triaxial	Yes**		

\*Coaxial sample holders require one FT-BNC or FT-TRIAX feedthrough as listed below

\*\*Triaxial sample holders require one FT-TRIAX feedthrough as listed below

\*\*\*Cannot be used in combination with the 360° sample stage rotation option

Part number	Description
FT-BNC	Coaxial feedthrough and coaxial cable, installed and wired
FT-TRIAX	Triaxial feedthrough and coaxial cable, installed and wired

### Site Requirements

**Power** — Standard single-phase (20 A maximum) power is required for instrumentation and optional pumps. The magnet power supply and optional recirculation chiller require 3-phase (21 A maximum) power.

**Water** — Electromagnet requires one supply and one return line for cooling with 7.6 L/min (2 gal/min) at a pressure drop of 30 psi. Magnet power supply requires 5.7 L/min (1.5 gal/min) at a pressure drop of 1.5 psi. Water temperature should be 15 °C to 25 °C.

## Equipment Options

Part Number	Description
PS-FOA	Optical fiber assembly. Transmit or receive light or IR/UV radiation. Fiber optic terminated with SMA connector or compression feedthrough. (Fiber optic wire not included. Fiber optic and probe cannot be used simultaneously.)
PS-Z16	16:1 zoom microscope upgrade; provides 8 µm resolution NOTE: consult Lake Shore for field upgrade
PS-LT	3.2 K base temperature stage pumping option—includes system modifications, stainless steel flex hose for pump connection, and Varian DS402 rotary pump (12.3 CFM at 60 Hz with oil mist eliminator on output [with oil return], inlet oil demister); 50 Hz operation may increase base temperature
PS-PLVI-25	Pump-line vibration isolator—includes NW 25 fittings, 1 m stainless steel bellows, clamps, and rings (placed between pump cart and system); requires one bag of cement. NOTE: for use with PT-LT

## Equipment Options, continued

<b>PS-360-EMPX</b>	360° sample stage rotation NOTE: limits low end temperature to 8 K; $\Delta T$ between the sample and the sample stage temperature sensor is $\leq \pm 2$ K over the full temperature range; recommended when taking measurements using microwave probes; consult Lake Shore for information on upgrading existing probe station
<b>PS-DPC</b>	Automatic Dewar (gas) pressure controller, regulates Dewar liquid flow
<b>PS-V81DPC</b>	Compact turbo pumping system—includes Varian V-81 turbo pump with oil free dry scroll backing pump, vacuum gauge for high vacuum, controller, and adaptors (specify 120 V/60 Hz or 220 V/50 Hz)
<b>PS-PLVI-40</b>	Pump-line vibration isolator—includes NW 40 fittings, 1 m stainless steel bellows, clamps, and rings (placed between pump cart and system); requires one bag of cement. NOTE: for use with PS-V81DPC
<b>PS-LN2</b>	Nitrogen Dewar with stainless fittings, gauges, and adaptors; allows LN <sub>2</sub> use with the LHe transfer line
<b>PS-OAC</b>	Oil-less compressor for PS-TMC (only available in 120 V)
<b>PA-SEN</b>	Probe arm modification with temperature sensor installed and wired to a 6-pin feedthrough
<b>PS-PAB-10</b>	Probe arm and base
<b>CS-5</b>	Calibration substrate for GSG probes—pad size: 50 $\mu\text{m}^2$ ; calibration type: SOLT (short-open-load-through), LRL (line-reflective-line), LRM (line-reflective-match); pitch range: 75 to 250 $\mu\text{m}$
<b>CS-15</b>	Calibration substrate for GSG probes—pad size: 25 $\mu\text{m}^2$ ; calibration type: SOLT (short-open-load-through), LRL (line-reflective-line), LRM (line-reflective-match); pitch range: 40 to 150 $\mu\text{m}$
<b>RC-EM7-208230-60-CE</b>	Recirculating chiller; P-2, 208 to 230 V, 60 Hz, 15 A, 5000 W, 15.1 L/min (4 gal/min), 60 psi, TF-5000
<b>RC-EM7-400-50-CE</b>	Recirculating chiller; P-2, 380 to 400 V, 50 Hz, 3-phase, 20 A, 6425 W, 12.5 L/min (3.3 gal/min), 60 psi, TF-7500
<b>RC-EM7-230-50-CE</b>	Recirculating chiller; P-2, 200 V, 50 Hz, 15 A, 4400 W, 12.5 L/min (3.3 gal/min), 60 psi, TF-5000

**Lake Shore Cryotronics is a leading supplier of cryogenic, superconducting magnet-based, electromagnet-based, high vacuum, and load-lock probe stations. We offer a full line of standard probe stations to meet your research requirements.**

*Our standard line of probe stations includes:*

### Model TTPX

- Up to six micro-manipulated probe stages
- Temperature range capabilities from 3.2 K to 675 K
- Up to 51 mm (2 in) diameter wafer capabilities
- Optical access through sample stage

### Model CPX

- Up to six micro-manipulated probe stages
- Temperature range capabilities from 1.5 K to 675 K
- Up to 51 mm (2 in) diameter wafer capabilities
- $\pm 5^\circ$  sample stage rotation
- Load-lock and high vacuum options

### Model FWPX

- Up to six micro-manipulated probe stages
- Temperature range capabilities from 3.5 K to 475 K
- $\pm 5^\circ$  sample stage rotation
- Up to 102 mm (4 in) diameter wafer capabilities

### Model CPX-HF

- 1 T horizontal field split pair superconducting magnet
- Up to four micro-manipulated probe stages
- Temperature range capabilities from 2 K to 400 K with field on or off
- $\pm 5^\circ$  sample stage rotation
- Up to 25 mm (1 in) diameter wafer capabilities
- High vacuum option

### Model CPX-VF

- 2.5 T vertical field solenoid superconducting magnet
- Up to six micro-manipulated probe stages
- Temperature range capabilities from 2 K to 400 K with field on or off
- $\pm 5^\circ$  sample stage rotation
- Up to 51 mm (2 in) diameter wafer capabilities
- High vacuum option

### Model EMPX-HF

- 0.55 T horizontal (in-plane field) electromagnet
- Up to four micro-manipulated probe stages
- Temperature range capabilities from 3.2 K to 400 K
- Up to 25 mm (1 in) diameter wafer capabilities
- 360° sample stage rotation option

### Model CRX-4K

- Low vibration, cryogen-free closed cycle refrigerator
- Up to six micro-manipulated probe stages
- Temperature range capabilities from 4.5 K to 475 K
- Up to 51 mm (2 in) diameter wafer capabilities