



Model 211 Temperature Controller



Model 211 features

- Operates down to 1.2 K with appropriate sensor
- One sensor input
- Supports diode and RTD sensors
- 0 V to 10 V or 4 mA to 20 mA output
- Large 5-digit LED display
- RS-232C serial interface and alarm relays
- CE certification
- Full 3 year standard warranty





Introduction

The Lake Shore single-channel Model 211 temperature monitor provides the accuracy, resolution, and interface features of a benchtop temperature monitor in an easy to use, easily integrated, compact instrument. With appropriate sensors, it measures from 1.2 K to 873 K, including temperatures in high vacuum and magnetic fields. Alarms, relays, user-configurable analog voltage or current output, and a serial interface are standard features on the Model 211. It is a good choice for liquefied gas storage and monitoring, cryopump control, cryo-cooler, and materials science applications, and when you need greater accuracy than thermocouples allow.

Sensor input reading capability

The Model 211 temperature monitor supports diode temperature sensors and resistance temperature detectors (RTDs). It can be configured for the type of sensor in use from the instrument front panel. Ensuring high accuracy and 5-digit measurement resolution are 4-lead differential measurement and 24-bit analog-to-digital conversion.

The Model 211 converts voltage or resistance to temperature units based on temperature response curve data for the sensor in use. Standard temperature response curves for silicon diodes and platinum RTDs are included in instrument firmware. It also provides non-volatile memory for one 200-point temperature response curve, which can be entered via the serial interface.

Interface

With an RS-232C serial interface and other interface features, the Model 211 is valuable as a stand-alone monitor and is easily integrated into other systems. Setup and every instrument function can be performed via serial interface or the front panel. Temperature data can be read up to seven times per second over computer interface; the display is updated twice each second. High and low alarms can be used in latching mode for error limit detection and in non-latching mode in conjunction with relays to perform simple on-off control functions. The analog output can be configured for either 0 to 10 V or 4 to 20 mA output.

Sensor Selection

Sensor temperature range (sensors sold separately)

| | | Model | Useful range | Magnetic field use |
|--|-----------------------|------------|-------------------------------|---|
| Diodes | Silicon diode | DT-670-SD | 1.4 K to 500 K | $T \geq 60 \text{ K} \ \& \ B \leq 3 \text{ T}$ |
| | Silicon diode | DT-670E-BR | 30 K to 500 K | $T \geq 60 \text{ K} \ \& \ B \leq 3 \text{ T}$ |
| | Silicon diode | DT-414 | 1.4 K to 375 K | $T \geq 60 \text{ K} \ \& \ B \leq 3 \text{ T}$ |
| | Silicon diode | DT-421 | 1.4 K to 325 K | $T \geq 60 \text{ K} \ \& \ B \leq 3 \text{ T}$ |
| | Silicon diode | DT-470-SD | 1.4 K to 500 K | $T \geq 60 \text{ K} \ \& \ B \leq 3 \text{ T}$ |
| | Silicon diode | DT-471-SD | 10 K to 500 K | $T \geq 60 \text{ K} \ \& \ B \leq 3 \text{ T}$ |
| | GaAlAs diode | TG-120-P | 1.4 K to 325 K | $T > 4.2 \text{ K} \ \& \ B \leq 5 \text{ T}$ |
| | GaAlAs diode | TG-120-PL | 1.4 K to 325 K | $T > 4.2 \text{ K} \ \& \ B \leq 5 \text{ T}$ |
| | GaAlAs diode | TG-120-SD | 1.4 K to 500 K | $T > 4.2 \text{ K} \ \& \ B \leq 5 \text{ T}$ |
| Positive temperature coefficient RTDs | 100 Ω platinum | PT-102/3 | 14 K to 873 K | $T > 40 \text{ K} \ \& \ B \leq 2.5 \text{ T}$ |
| | 100 Ω platinum | PT-111 | 14 K to 673 K | $T > 40 \text{ K} \ \& \ B \leq 2.5 \text{ T}$ |
| | Rhodium-iron | RF-800-4 | 1.4 K to 500 K | $T > 77 \text{ K} \ \& \ B \leq 8 \text{ T}$ |
| | Rhodium-iron | RF-100T/U | 1.4 K to 325 K | $T > 77 \text{ K} \ \& \ B \leq 8 \text{ T}$ |
| Negative temperature coefficient RTDs ¹ | Cernox™ | CX-1010 | 2 K to 325 K ⁴ | $T > 2 \text{ K} \ \& \ B \leq 19 \text{ T}$ |
| | Cernox™ | CX-1030-HT | 3.5 K to 420 K ^{2,5} | $T > 2 \text{ K} \ \& \ B \leq 19 \text{ T}$ |
| | Cernox™ | CX-1050-HT | 4 K to 420 K ^{2,5} | $T > 2 \text{ K} \ \& \ B \leq 19 \text{ T}$ |
| | Cernox™ | CX-1070-HT | 15 K to 420 K ² | $T > 2 \text{ K} \ \& \ B \leq 19 \text{ T}$ |
| | Cernox™ | CX-1080-HT | 50 K to 420 K ² | $T > 2 \text{ K} \ \& \ B \leq 19 \text{ T}$ |
| | Germanium | GR-300-AA | 1.2 K to 100 K ³ | Not recommended |
| | Germanium | GR-1400-AA | 4 K to 100 K ³ | Not recommended |
| | Rox™ | RX-102A | 1.4 K to 40 K ⁴ | $T > 2 \text{ K} \ \& \ B \leq 10 \text{ T}$ |

¹ Single excitation current may limit the low temperature range of NTC resistors

² Non-HT version maximum temperature: 325 K

³ Low temperature limited by input resistance range

⁴ Low temperature specified with self-heating error: $\leq 5 \text{ mK}$


⁵ Low temperature specified with self-heating error: $\leq 12 \text{ mK}$

Silicon diodes are the best choice for general cryogenic use from 1.4 K to above room temperature. Diodes are economical to use because they follow a standard curve and are interchangeable in many applications. They are not suitable for use in ionizing radiation or magnetic fields.

Cernox™ thin-film RTDs offer high sensitivity and low magnetic field-induced errors over the 2 K to 420 K temperature range. Cernox sensors require calibration.

Platinum RTDs offer high uniform sensitivity from 30 K to over 800 K. With excellent reproducibility, they are useful as thermometry standards. They follow a standard curve above 70 K and are interchangeable in many applications.

Model 211 rear panel



1 Power input connector

2 Serial (RS-232C) I/O (DTE)

3 Analog output



Display

The Model 211 has a 6-digit LED display with measurements available in temperature units K, °C, °F, or sensor units V or Ω .

Specifications

Sensor input configuration

| Diode/RTD | |
|-------------------|--|
| Measurement type | 4-lead differential |
| Excitation | 8 constant current sources |
| Supported sensors | Diodes: silicon, GaAlAs RTDs: 100 Ω platinum, 1000 Ω platinum, germanium, carbon-glass, Cernox™, and Rox™ |
| Standard curves | DT-470, DT-670, CTI-C, PT-100, and PT-1000 |
| Input connector | Shared 25-pin D-sub |

Thermometry

Number of inputs 1

Input configuration Input can be configured from the front panel to accept any of the supported input types

Isolation Measurement is not isolated from chassis ground

A/D resolution 24-bit

Input accuracy Sensor dependent—refer to Input Specifications table

Measurement resolution Sensor dependent—refer to Input Specifications table

Maximum update rate 7 rdg/s

User curve One 200-point CalCurve™ or user curve in non-volatile memory

Front panel

Display 5-digit LED

Number of reading displays 1

Display units K, °C, °F, V, and Ω

Reading source Temperature and sensor units

Display update rate 2 rdg/s

Temp display resolution 0.001° from 0° to 99.999°, 0.01° from 100° to 999.99°, 0.1° above 1000°

Sensor units display resolution Sensor dependent to 5 digits

Display annunciators K, °C, °F, and V/ Ω

Keypad 4 full travel keys, numeric and specific functions

Front panel features Display brightness control, keypad lock-out

Typical sensor performance—see Appendix F for sample calculations of typical sensor performance

| | Example Lake Shore sensor | Temperature | Nominal resistance/voltage | Typical sensor sensitivity ⁶ | Measurement resolution: temperature equivalents | Electronic accuracy: temperature equivalents | Temperature accuracy including electronic accuracy, CalCurve™, and calibrated sensor |
|--|--|-------------|----------------------------|---|---|--|--|
| Silicon diode | DT-670-SD with 1.4H calibration | 1.4 K | 1.644 V | -12.49 mV/K | 1.6 mK | ±29 mK | ±41 mK |
| | | 77 K | 1.028 V | -1.73 mV/K | 11.6 mK | ±175 mK | ±197 mK |
| | | 300 K | 0.5597 V | -2.3 mV/K | 8.7 mK | ±111 mK | ±143 mK |
| | | 500 K | 0.0907 V | -2.12 mV/K | 9.4 mK | ±99 mK | ±149 mK |
| Silicon diode | DT-470-SD-13 with 1.4H calibration | 1.4 K | 1.6981 V | -13.1 mV/K | 1.5 mK | ±28 mK | ±40 mK |
| | | 77 K | 1.0203 V | -1.92 mV/K | 10.5 mK | ±157 mK | ±179 mK |
| | | 300 K | 0.5189 V | -2.4 mV/K | 8.4 mK | ±105 mK | ±137 mK |
| | | 475 K | 0.0906 V | -2.22 mV/K | 9.1 mK | ±94 mK | ±144 mK |
| GaAlAs diode | TG-120-SD with 1.4H calibration | 1.4 K | 5.391 V | -97.5 mV/K | 0.2 mK | ±15 mK | ±27 mK |
| | | 77 K | 1.422 V | -1.24 mV/K | 16.2 mK | ±512 mK | ±534 mK |
| | | 300 K | 0.8978 V | -2.85 mV/K | 7 mK | ±186 mK | ±218 mK |
| | | 475 K | 0.3778 V | -3.15 mV/K | 6.4 mK | ±135 mK | ±185 mK |
| 100 Ω platinum RTD 500 Ω full scale | PT-103 with 1.4J calibration | 30 K | 3.66 Ω | 0.19 Ω /K | 10.5 mK | ±320 mK | ±330 mK |
| | | 77 K | 20.38 Ω | 0.42 Ω /K | 4.8 mK | ±153 mK | ±165 mK |
| | | 300 K | 110.35 Ω | 0.39 Ω /K | 5.2 mK | ±210 mK | ±232 mK |
| | | 500 K | 185.668 Ω | 0.378 Ω /K | 5.3 mK | ±257 mK | ±303 mK |
| Cernox™ | CX-1050-SD-HT ⁷ with 4M calibration | 4.2 K | 3507.2 Ω | -1120.8 Ω /K | 45 μ K | ±2.0 mK | ±7.0 mK |
| | | 77 K | 205.67 Ω | -2.4116 Ω /K | 20.8 mK | ±366 mK | ±382 mK |
| | | 300 K | 59.467 Ω | -0.1727 Ω /K | 290 mK | ±4.8 K | ±4.8 K |
| | | 420 K | 45.03 Ω | -0.0829 Ω /K | 604 mK | ±9.9 K | ±9.9 K |
| Germanium | GR-300-AA with 0.3D calibration | 1.2 K | 600 Ω | -987 Ω /K | 51 μ K | ±0.6 mK | ±5.3 mK |
| | | 1.4 K | 449 Ω | -581 Ω /K | 86 μ K | ±1 mK | ±5 mK |
| | | 4.2 K | 94 Ω | -27 Ω /K | 1.9 mK | ±16 mK | ±20 mK |
| | | 100 K | 3 Ω | -0.024 Ω /K | 2.10 K | ±2.5 K | ±2.5 K |
| Germanium | GR-1400-AA with 1.4D calibration | 4 K | 1873 Ω | -1008 Ω /K | 50 μ K | ±1.1 mK | ±5.1 mK |
| | | 4.2 K | 1689 Ω | -862 Ω /K | 58 μ K | ±1.2 mK | ±5.2 mK |
| | | 10 K | 253 Ω | -62 Ω /K | 807 μ K | ±1.8 mK | ±6.3 mK |
| | | 100 K | 3 Ω | -0.021 Ω /K | 2.40 K | ±2.9 K | ±2.9 K |
| Carbon-glass (no longer available) | CGR-1-2000 with 4L calibration | 4.2 K | 2260 Ω | -2060 Ω /K | 25 μ K | ±0.6 mK | ±4.6 mK |
| | | 77 K | 21.65 Ω | -0.157 Ω /K | 319 mK | ±410 mK | ±435 mK |
| | | 300 K | 11.99 Ω | -0.015 Ω /K | 3.33 K | ±4.2 K | ±4.2 K |

⁶ Typical sensor sensitivities were taken from representative calibrations for the sensor listed

⁷ Non-HT version maximum temperature: 325 K



Input specifications

| Sensor type | Sensor temperature coefficient | Input range | Excitation current | Display resolution | Measurement resolution | Electronic accuracy | Instrument temperature coefficient |
|--|--------------------------------|-----------------------------|-------------------------------------|--------------------|------------------------|--|---|
| Silicon diode | negative | 0 V to 2.5 V | 10 μ A \pm 0.05% ⁸ | 100 μ V | 20 μ V | \pm 200 μ V \pm 0.01% of rdg | \pm 10 μ V \pm 5 PPM of rdg/ $^{\circ}$ C |
| GaAlAs diode | negative | 0 V to 7.5 V | 10 μ A \pm 0.05% ⁸ | 100 μ V | 20 μ V | \pm 350 μ V \pm 0.02% of rdg | \pm 20 μ V \pm 5 PPM of rdg/ $^{\circ}$ C |
| 100 Ω platinum RTD, 250 Ω full scale | positive | 0 Ω to 250 Ω | 1 mA \pm 0.3% ⁹ | 10 m Ω | 2 m Ω | \pm 0.06 Ω \pm 0.02% of rdg | \pm 0.2 m Ω \pm 5 PPM of rdg/ $^{\circ}$ C |
| 100 Ω platinum RTD, 500 Ω full scale | positive | 0 Ω to 500 Ω | 1 mA \pm 0.3% ⁹ | 10 m Ω | 2 m Ω | \pm 0.06 Ω \pm 0.02% of rdg | \pm 0.2 m Ω \pm 5 PPM of rdg/ $^{\circ}$ C |
| 1000 Ω platinum RTD | positive | 0 Ω to 5000 Ω | 1 mA \pm 0.3% ⁹ | 100 m Ω | 20 m Ω | \pm 0.4 Ω \pm 0.04% of rdg | \pm 2.0 m Ω \pm 5 PPM of rdg/ $^{\circ}$ C |
| Cernox™ RTD | negative | 0 Ω to 7500 Ω | 10 μ A \pm 0.05% ⁸ | 100 m Ω | 50 m Ω | \pm 0.8 Ω \pm 0.04% of rdg | \pm 20 m Ω \pm 15 PPM of rdg/ $^{\circ}$ C |

⁸ Current source error has negligible effect on measurement accuracy

⁹ Current source error is removed during calibration

Interface

Serial interface

| | |
|--------------------------|---------------|
| Electrical format | RS-232C |
| Max baud rate | 9600 baud |
| Connector | 9-pin D-sub |
| Reading rate | Up to 7 rdg/s |

Alarms

| | |
|--------------------|--|
| Number | 2, high and low |
| Data source | Temperature |
| Settings | High setpoint, Low setpoint, Dead band, Latching or Non-latching |
| Actuators | Display message, relays |

Relays

| | |
|-----------------------|--|
| Number | 2 |
| Contacts | Normally Open (NO), Normally Closed (NC), and Common (C) |
| Contact rating | 30 VDC at 1 A |
| Operation | Activate relays on high or low input alarm or manual |
| Connector | Shared 25-pin D-sub |

Analog output

| | |
|--------------------|--|
| Isolation | Output is not isolated from chassis ground |
| Update rate | 7 readings per s |
| Data source | Temperature |

| | Voltage | Current |
|----------------------------|---------------|--|
| Range | 0 V to 10 V | 4 mA to 20 mA |
| Accuracy | \pm 1.25 mV | \pm 5.0 μ A |
| Resolution | 0.3 mV | 0.6 μ A |
| Min load resistance | 500 Ω | NA |
| Compliance voltage | NA | 10 V |
| Load regulation | NA | \pm 0.02% of reading 0 to 500 Ω |

| | Temperature | Sensor units (fixed by type) |
|---------|---------------|---|
| Scales: | 0 K to 20 K | Diodes: 1 V = 1 V |
| | 0 K to 100 K | 100 Ω platinum: 1 V = 100 Ω |
| | 0 K to 200 K | 1000 Ω platinum: 1 V = 1000 Ω |
| | 0 K to 325 K | NTC resistor: 1 V = 1000 Ω |
| | 0 K to 475 K | |
| | 0 K to 1000 K | |

Settings Voltage or current, scale

Connector Shared 25-pin D-sub

General

Ambient temperature 15 $^{\circ}$ C to 35 $^{\circ}$ C at rated accuracy, 10 $^{\circ}$ C to 40 $^{\circ}$ C at reduced accuracy

Power requirements Regulated +5 VDC at 400 mA

Size 96 mm W \times 48 mm H \times 166 mm D (3.8 in \times 1.9 in \times 6.5 in)

Mounting Panel mount into 91 mm W \times 44 mm H (3.6 in \times 1.7 in) cutout

Weight 0.45 kg (1 lb)

Approvals CE mark, RoHS



2111 Single 1/4 DIN panel-mount adapter, 105 mm W \times 132 mm H (4.1 in \times 5.2 in)



2112 Dual 1/4 DIN panel-mount adapter, 105 mm W \times 132 mm H (4.1 in \times 5.2 in)

Power supply (109-132)


Comes standard with interchangeable input plugs

Power requirements 100 to 240 VAC, 50 or 60 Hz, 0.3 A max

Output +5 V at 1.2 A

Size 40.5 mm W \times 30 mm H \times 64 mm D (1.6 in \times 1.2 in \times 2.5 in)

Weight 0.15 kg (0.33 lb)



Ordering information

| Part number | Description |
|-------------|--|
| 211S | Model 211 single channel temperature monitor—includes 100 to 240 V, 6 W universal power supply with interchangeable input plugs (109-132), one DB-25 sensor input mating connector (G-106-253), one sensor input mating connector shell (G-106-264), a calibration certificate and a user's manual |
| 211N | Model 211N with all accessories except the power supply |

Accessories

| | |
|---------------------|--|
| 109-132 | 100-240 VAC power supply with interchangeable plugs for US, UK, Europe, Australia, and China application |
| 2111 | Single 1/4 DIN panel-mount adapter |
| 2112 | Dual 1/4 DIN panel-mount adapter |
| 8000 | CalCurve™, CD-ROM (included with calibrated sensor) |
| G-106-253 | DB-25 plug, qty 1 |
| G-106-264 | DB-25 hood, qty 1 |
| CAL-211-CERT | Instrument recalibration with certificate |
| CAL-211-DATA | Instrument recalibration with certificate and data |
| 119-043 | Model 211 temperature monitor manual |

All specifications are subject to change without notice

