



CryoComplete™

Spanning the cryogenic ecosystem

77 K to 500 K

Everything you need
to start making
temperature-dependent,
low-level electrical
measurements



environment by  JANIS

Focus your research with CryoComplete

With CryoComplete™, you can start making cryogenic electrical measurements as soon as it lands in your lab. From the simple-to-use, pour-fill LN₂ Dewar to the prewritten I-V (resistance) measurement routines, CryoComplete produces results right out of the box.

While easy to use, the system's performance doesn't disappoint. Its industry-leading measurement electronics promote low-level DC measurements and three full channels of lock-in AC capability—the keys to unlocking difficult measurements. Best of all, our cryogenic experts have designed CryoComplete from top to bottom, using cryogenic best practices, to deliver end-to-end system specifications.

PC with MeasureLINK™

A PC with MeasureLINK provides the user interface to control your cryogenic system. MeasureLINK enables a wide range of capabilities, including data charting, instrument control, and system monitoring with a cryostat-specific process view.

LN₂ cryostat

The Environment by Janis VPF-100 sample in vacuum cryostat provides a variable-temperature sample environment from 77 K to 500 K. The pour-fill design allows quick and easy LN₂ refills.



Source + measure + lock-in

Run ultra-low-noise AC/DC measurements with the MeasureReady™ M81-SSM synchronous source measure system. In addition to the M81-SSM-6 instrument, it includes a BCS-10 balanced current source module and a VM-10 DC/AC/lock-in voltmeter module with a combined noise performance (differential) of 4.1 nV/√Hz.

Temperature control


Control temperature within 50 mK with a Lake Shore Model 335 temperature controller, a Lake Shore precision-calibrated silicon diode, and a pre-wired heater. Advanced PID autotuning, pre-programmed sensor calibration, and default cryostat tuning enable fast setup and operation.

Applications and capabilities


From setup to measurement, CryoComplete enhances your cryogenic experimentation.

- Complete measurement system
- Optimized signal path
- Quick lead times

Measurement benefits



Common measurements



Thermal transport	1D materials, thermoelectric materials	✓	✓					✓				✓
Materials research	Nanodevices, superconducting devices, nonlinear devices			✓				✓	✓	✓	✓	
Materials development	Linear systems, sensors			✓						✓	✓	

Standard system capabilities

VPF-100 cryostat/335 temperature controller/calibrated silicon diode

- Operating temperature range:** 77 K to 500 K
Cryogen: Liquid nitrogen
Sample environment: Sample in vacuum
Temperature stability: 50 mK
Pour-fill reservoir capacity: 1.2 L LN₂
Cooldown time: 15 min to 77 K
Working time: 8 h
Optical ports: 4 quartz windows
Electrical sample mount: Pre-wired mounting plate with 8 contact pins

Resistance/I-V measurements

M81-SSM-6 with balanced current source and voltmeter modules


- Measurements:** 100 $\mu\Omega$ to 1 G Ω *
Source modes: DC, sine, triangle, square
Source ranges: 1 pA to 100 mA
Source frequency: 100 μ Hz to 100 kHz (square <5 kHz)

**Upper impedance range limited to DC*

- Measurement limit:** 10 V maximum
Input impedance: >10 G Ω (differential)

- Leakage current at sample:** 50 pA at 10 V for coaxial or 50 fA at 10 V for guarded triaxial
Voltage noise at sample: <5 nV/ $\sqrt{\text{Hz}}$ at 83 Hz
Measure noise at sample (1/f): <100 nV


BCS-10 versus VM-10, 10 m Ω resistor, 4-probe, 2TX and 2CXLIA at 83 Hz, FIR = 3, τ = 200 ms





Ordering information

CryoComplete-LN2-V
CryoComplete 77 K to 500 K cryogenic characterization system


Easily control and monitor your system with **MeasureLINK™** software



◀ **Process view** shows a representation of the cryostat internals with the appropriate temperatures highlighted for a better understanding of internal temperature variations (shown is an internal view of a VPF-100 application)




Create multiple measurement configurations



Monitor pane

The monitor pane allows easy access to monitor temperature and change control setpoints



Create nested, multi-level measurement loop sequences with drag-and-drop controls, and coordinate the cryostat environment with electrical source sweeps and multi-channel data collection