



Faster Hall analysis in a Quantum Design PPMS®
 Measure**L&NK**™ M91 FastHall™ controller

The Lake Shore Cryotronics M91 instrument delivers significantly higher levels of precision, speed, and convenience as compared to traditional Hall effect measurement solutions. Combining all the necessary Hall measurement functions into a single instrument, the M91 automatically executes measurements and calculates the final Hall and mobility parameters. The instrument eliminates manual trial-and-error steps and ensures that measurements are always made under optimal conditions for the sample.

The instrument also provides better measurements faster, especially when characterizing low-mobility materials. Most commonly measured materials can be analyzed in a few seconds. This is because the instrument's patented FastHall technology eliminates the need to reverse the magnetic field during the measurement. This is particularly beneficial when using the M91 with the superconducting magnet of a PPMS, allowing for much faster analysis of low-mobility material samples in van der Pauw (vdP) wired geometries.

Key features

No need for field reversal during the measurement due to the FastHall technology

Automatically selects optimal excitation and measurement levels, checks sample contact quality, executes measurement steps with standard protocols, and provides complete Hall analysis

Calculates derived parameters for vdP and Hall bar samples

FastHall technology extends mobility range down to 0.001 cm²/(V·s) when measuring vdP samples

Touchscreen UI displays measurement process steps as they execute in real time



The M91 automatically checks and graphically displays sample contact quality



Electrical measurement specifications

The M91 FastHall Measurement Controller integrates all the required source measure, and signal switching capabilities for these system measurement specifications:

Resistance (R)

Range: 10 mΩ to 10 MΩ source current

Monitor temperature and change setpoints with the monitor pane

Mobility (μ)

Mobility Range: 10⁶ cm²/(V·s) to 0.001 cm²/(V·s)

Source parameters


Current Source Range: 1 μA to 100 mA (lowest usable current: 10 nA)

Current Measurement Range: 100 mA to 10 nA

Compliance Voltage: 10 V, maximum


Voltage Measurement Range: 1 mV to 10 V

Compliance Current: 100 mA, maximum



For more information, visit www.lakeshore.com

Simple setup and integration



M91 to PPMS connections are made via triaxial cables available from Lake Shore Cryotronics. Samples are wired to the PPMS sample puck for either vdP (4-wire) or Hall bar (6-wire) geometries. During general operation, the sample connector is connected to the PPMS electronics, and the system's MultiVu™ software controls the field and temperature of the sample space while coordinating the measurement electronics.

Electrical measurement specifications

The M91's MeasureLINK™ software for Hall measurement control integrates easily with the PPMS MultiVu application software. MeasureLINK can be installed on the same PC with the MultiVu software or on a separate PC that is on the same network as the MultiVu PC.

Once installed, MeasureLINK:


Provides a simple way to start and step through your measurement sequences, as well as chart, log, and organize the results

Includes scripts for running Hall measurements and reporting the results

Enables automated control of field and management of sample temperature

Generates detailed reports including all the supporting intermediate data so you can readily confirm the integrity of the final results

Allows for customization of measurement sequences for specific Hall research requirements (optional upgrade)



Field-dependent transverse and longitudinal transport measurements for a GaAs 2-D electron gas system at 2 K with 1 μ A sourced excitation current in the van der Pauw geometry. Plateaux in the transverse channel demonstrate the integer quantum Hall effect and correspond to where the Fermi level falls in an area of localized states between neighboring Landau levels.

Sample provided by Dr. M. Pendharkar, Chris Palmström Group, University of California Santa Barbara.



www.qdusa.com • APS Booth #214

