


## TECHNICAL NOTE

# Real-Time FORC (RTForc™) Software for the 8600 Series VSM

B. C. Dodrill, H. S. Reichard, and T. Shimizu, Lake Shore Cryotronics



### Key Features

- Fully automated FORC data acquisition using the 8600 VSM software
- FORC distributions are calculated and displayed in real-time, significantly reducing the time required to collect and analyze FORC data
- Users can change between the ( $H_a$ ,  $H_b$ ) and ( $H_c$ ,  $H_u$ ) coordinate systems, select smoothing factors (SF) and number of contours to be displayed, export the FORC diagram image, etc.
- Output data is compatible with FORCinel post-processing software

## Background

First-order-reversal-curve (FORC) measurements and analysis provide information regarding magnetic reversal mechanisms in magnetic materials that cannot be obtained from major hysteresis loop measurements alone. It has been extensively used by earth and planetary scientists studying the magnetic properties of natural samples (rocks, soils, sediments, etc.) because FORC can distinguish between single-domain (SD), multi-domain (MD), and pseudo single-domain (PSD) behavior, and because it can discriminate between different magnetic mineral species<sup>1, 2</sup>. It has proven to be useful in better understanding the nature of magnetization reversal and interactions in magnetic nanowires<sup>3-7</sup>, nanomagnet arrays<sup>8-11</sup>, thin film magnetic recording media<sup>12-14</sup> and thin film magnetic multilayers<sup>15-17</sup>, nanostructured permanent magnet materials<sup>18,19</sup>, soft magnetic bilayers<sup>20</sup> and magneto-caloric effect (MCE) materials<sup>21</sup>. It has also been used to differentiate between phases in multiphase magnetic materials because it is very difficult to unravel the complex magnetic signatures of such materials from a hysteresis loop measurement alone<sup>22-24</sup>.

A FORC is measured by saturating a sample in a field  $H_{sat}$ , decreasing the field to a reversal field  $H_a$ , then measuring moment versus field  $H_b$  as the field is swept back to  $H_{sat}$ . This process is repeated for many values of  $H_a$ , yielding a series of FORCs as shown in figure 1 for a magnetic tape. The FORC distribution  $\rho(H_a, H_b)$  is the mixed second derivative:

$$\rho(H_a, H_b) = -(1/2)\partial^2 M(H_a, H_b)/\partial H_a \partial H_b$$

A FORC diagram is a 2D or 3D contour plot of  $\rho(H_a, H_b)$ . It is common to change the coordinates from  $(H_a, H_b)$  to:

$$H_c = (H_b - H_a)/2, H_u = (H_b + H_a)/2$$

$H_u$  represents the distribution of interaction or reversal fields, and  $H_c$  represents the distribution of switching or coercive fields. The 2D FORC diagram for the magnetic tape is shown in figure 2.

A FORC diagram not only provides information regarding the distribution of interaction and switching fields, but also serves as a “fingerprint” that gives insight into the domain state and nature of interactions occurring in magnetic materials. In a FORC diagram entirely closed contours are usually associated with SD behavior, while open contours that diverge towards the  $H_u$  axis are associated with MD, and open and closed contours together are associated with PSD. The peak in the FORC distribution is usually centered at a switching field  $H_c$  that correlates with the coercivity as

determined from a hysteresis loop measurement. If the peak in the FORC distribution is centered at an interaction field  $H_u = 0$  this means that interactions are weak. Conversely, if the peak is shifted towards positive  $H_u$  they are strong. Multiple peaks in a FORC diagram mean there are multiple magnetic phases in a material. And the very shape of the FORC distribution provides insight into the nature of interactions (dipolar, exchange) that are occurring in a magnetic material.

## The use of RTForc™ for FORC measurements

Conducting FORC measurements and analysis is usually an iterative process whereby one defines the FORC data acquisition parameters (figure 3 shows an example FORC Measurement setup), conducts the measurement, and then processes the data to generate a FORC diagram using 3<sup>rd</sup> party software (e.g., FORCinel<sup>25</sup>, VARIFORC<sup>26</sup>). After the first FORC measurement, adjustments are often made to the FORC acquisition parameters and the sample is re-measured to better optimize the FORC diagram, e.g., ( $H_c$ ,  $H_u$ ) range, resolution, etc. Since FORC measurements can take tens of minutes to several hours or longer (depending on the density of data being collected, and the sample being measured), the time required to produce the desired FORC diagram can be quite long. Real-time FORC generates the FORC diagram as data is being collected, enabling one to easily see if FORC parameters were defined correctly, and thus significantly reduces the amount of time required to optimally collect and analyze FORC data.




Figure 1: Measured FORCs for a magnetic tape.

## Example RTForc diagrams:




Figure 2: RTForc diagram for the magnetic tape. The peak is centered at  $H_c$  which correlates with the coercivity of the sample as determined from a hysteresis loop measurement. The peak is shifted towards negative interaction fields  $H_u$  and the distribution has a “boomerang” shape. These are features that are usually associated with exchange interactions.




Figure 4: RTForc diagram for a marine sediment sample. The FORC diagram consists of both closed and open contours, the latter diverging towards the  $H_u$  axis. These features are usually associated with pseudo-single domain (PSD) behavior.

The screenshot shows the FORC Measurement software interface. It is divided into several sections:
 

- Saturation field:** Saturation field: 10 kOe, Default FORC settings button.
- FORC Spacing:** Number of FORCs: 75, Pause at saturation field: 0.10 seconds, Pause at reversal fields: 1.00 seconds, Field step size: 40.540541 Oe, Averaging time: 0.10 seconds, Estimated execution time (hh:mm:ss): 00:15:26.
- FORC data acquisition range:** Classic mode selected. Minimum  $H_u$  field: -750 Oe, Maximum  $H_u$  field: 750 Oe, Maximum  $H_c$  field: 1.5 kOe, Pause at calibration field: 1.00 seconds. Extended range mode is unselected.
- FORC Analysis (Classic mode only):** Show FORC diagram on the fly (checked), Smoothing factor: 3, Rotate 45 degrees (checked), Truncate for a rectangular diagram (checked), Number of contours: 15, Enable re-gridding (checked), Correct for drift (checked).

Figure 3: Example 8600 VSM FORC Measurement setup.



Watch a video demonstrating an RTForc™ measurement




Figure 5: RTForc diagram for a  $BaFe_{12}O_{19}$  nanoparticles (NP) sample. There are two peaks corresponding to low (soft) and high (hard) coercivity components, and the region between the two peaks is related to the exchange coupling between the two phases<sup>24</sup>.




Figure 6: RTForc diagram for a permanent magnet sample. The FORC diagram has a “wishbone” like feature which is usually associated with long-range magneto-static (dipolar) interactions.

## References:

1. C. R. Pike, A. P. Roberts, and K. L. Verosub, Characterizing Interactions in Fine Magnetic Particle Systems Using First Order Reversal Curves, *J. Appl. Phys.*, 85, 6660, 1999.
2. A. P. Roberts, C. R. Pike, and K. L. Verosub, First-Order Reversal Curve Diagrams: A New Tool for Characterizing the Magnetic Properties of Natural Samples, *J. Geophys. Res.*, 105, 461, 2000.
3. A. Rotaru, J. Lim, D. Lenormand, A. Diaconu, J. Wiley, P. Postolache, A. Stancu, and L. Spinu, Interactions and Reversal Field Memory in Complex Magnetic Nanowire Arrays, *Phys. Rev. B*, 84, 13, 134431, 2011.
4. O. Trusca, D. Cimpoesu, J. Lim, X. Zhang, J. Wiley, A. Diaconu, I. Dumitru, A. Stancu, and L. Spinu, Interaction Effects in Ni Nanowire Arrays, *IEEE Trans. Mag.*, 44, 11, 2730, 2008.
5. A. Arefpour, M. Almasi-Kashi, A. Ramazani, and E. Golafshan, The Investigation of Perpendicular Anisotropy of Ternary Alloy Magnetic Nanowire Arrays Using First Order Reversal Curves, *J. Alloys and Comp.*, 583, 340, 2014.
6. B. C. Dodrill, and L. Spinu, First-Order-Reversal-Curve Analysis of Nanoscale Magnetic Materials, Technical Proceedings of the 2014 NSTI Nanotechnology Conference and Exposition, CRC Press, 2014.
7. A. Sharma, M. DiVito, D. Shore, A. Block, K. Pollock, P. Solheid, J. Feinberg, J. Modiano, C. Lam, A. Hubel, and B. Stadler, Alignment of Collagen Matrices Using Magnetic Nanowires and Magnetic Barcode Readout Using First Order Reversal Curves (FORC), *J. Mag. Mag. Mat.*, 459, 176, 2018.
8. F. Beron, L. Carignan, D. Menard, and A. Yelon, Extracting Individual Properties from Global Behavior: First Order Reversal Curve Method Applied to Magnetic Nanowire Arrays, *Electrodeposited Nanowires and Their Applications*, edited by N. Lupu, 228, INTECH, Croatia, 2010.
9. R. Dumas, C. Li, I. Roshchin, I. Schuller, and K. Liu, Magnetic Fingerprints of sub-100 nm Fe Nanodots, *Phys. Rev. B*, 75, 134405, 2007.
10. D. Gilbert, G. Zimanyi, R. Dumas, M. Winklhofer, A. Gomez, N. Eibagi, J. Vincent, and K. Liu, Quantitative Decoding of Interactions in Tunable Nanomagnet Arrays, *Sci. Reports*, 4, 4204, 2014.
11. J. Graffe, M. Weigand, C. Stahl, N. Trager, M. Kopp, G. Schutz, and E. Goering, Combined First Order Reversal Curve and X-ray Microscopy Investigation of Magnetization Reversal Mechanisms in Hexagonal Antidot Arrays, *Phys. Rev. B*, 93, 014406, 2016.
12. B. Valcu, D. Gilbert, and K. Liu, Fingerprinting Inhomogeneities in Recording Media Using First Order Reversal Curves, *IEEE Trans. Mag.*, 47, 2988, 2011.
13. A. Stancu and E. Maccim, Interaction Field Distribution in Longitudinal and Perpendicular Structured Particulate Media, *IEEE Trans. Mag.*, 42, 10, 3162, 2006.
14. M. Winklhofer, R. K. Dumas, and K. Liu, Identifying Reversible and Irreversible Magnetization Changes in Prototype Patterned Media Using First- and Second-Order Reversal Curves, *J. Appl. Phys.*, 103, 07C518, 2008.
15. R. Dumas, C. Li, L. Roshchin, I. Schuller, and K. Liu, Deconvoluting Reversal Modes in Exchange Biased Nanodots, *Phys. Rev. B*, 144410, 2012.
16. R. Gallardo, S. Khanai, J. Vargas, L. Spinu, C. Ross, and C. Garcia, Angular Dependent FORC and FMR of Exchange Biased NiFe Multilayer Films, *J. Phys. D: Appl. Phys.*, 50, 075002, 2017.
17. N. Siadou, M. Androutsopoulos, I. Panagiotopoulos, L. Stoleriu, A. Stancu, T. Bakas, and V. Alexandrakis, Magnetization Reversal in  $[\text{Ni}/\text{Pt}]_6/\text{Pt}(x)/[\text{Co}/\text{Pt}]_6$  Multilayers, *J. Mag. Mag. Mat.*, 323, 12, 1671, 2011.
18. T. Schrefl, T. Shoji, M. Winklhofer, H. Oezeit, M. Yano, and G. Zimanyi, First Order Reversal Curve Studies of Permanent Magnets, *J. Appl. Phys.*, 111, 07A728, 2012.
19. M. Pan, P. Shang, H. Ge, N. Yu, and Q. Wu, First Order Reversal Curve Analysis of Exchange Coupled SmCo/NdFeB Nanocomposite Alloys, *J. Mag. Mag. Mat.*, 361, 219, 2014.
20. M. Rivas, J. Garcia, I. Skorvanek, J. Marcin. P. Svec, and P. Gorria, Magnetostatic Interaction in Soft Magnetic Bilayer Ribbons Unambiguously Identified by First Order Reversal Curve Analysis, *Appl. Phys. Letts.*, 107, 132403, 2015.
21. V. Franco, F. Beron. K. Pirota, M. Knobel, and M. Willard, Characterization of Magnetic Interactions of Multiphase Magnetocaloric Materials using First Order Reversal Curve Analysis, *J. Appl. Phys.*, 117, 17C124, 2015.
22. B. C. Dodrill, First-Order-Reversal-Curve Analysis of Nanocomposite Permanent Magnets, Technical Proceedings of the 2015 TechConnect World Innovation Conference and Expo, CRC Press, 2015.
23. C. Carvallo, A. R. Muxworthy, and D. J. Dunlop, First-Order-Reversal-Curve (FORC) Diagrams of Magnetic Mixtures: Micromagnetic Models and Measurements, *Physics of the Earth and Planetary Interiors*, 154, 308, 2006.
24. Y. Cao, M. Ahmadzadeh, K. Xe, B. Dodrill, and J. McCloy, Multiphase Magnetic Systems: Measurement and Simulation, *J. Appl. Phys.*, 123(2), 023902, 2018.
25. R. J. Harrison and J. M. Feinberg, FORCinel: An Improved Algorithm for Calculating First-Order Reversal Curve Distributions Using Locally Weighted Regression Smoothing, *Geochemistry, Geophysics, Geosystems*, 9, 11, 2008.
26. R. Egli, VARIFORC: An Optimized Protocol for Calculating Non-regular First-Order Reversal Curve (FORC) Diagrams, *Global and Planetary Change*, 203, 110, 203, 2013.