




# 8600 Series VSM Measurement Results

B. C. Dodrill





## Typical Noise and Low Moment Measurement Results

### NOISE AND LOW MOMENT

#### Figure 1a and 1b


Noise at 100 ms/point (1a) and 10 s/point (1b) averaging. The observed noise is 119.5 nemu and 13 nemu RMS, respectively, consistent with our RMS noise specifications of 150 nemu and 15 nemu, respectively. Measurements are taken with an empty sample holder.




### NOISE AND LOW MOMENT

#### Figure 2

1 min 25 s at 100 ms/point hysteresis loop for a 20  $\mu$ emu CoPt bit-patterned magnetic media (thin film) sample. The peak-to-peak noise in the data is consistent with the RMS noise specification at 100 ms/point.






## NOISE AND LOW MOMENT

### Figure 3


13 min 30 s at 1 s/point hysteresis loop for the 20  $\mu\text{emu}$  CoPt bit-patterned magnetic media (thin film) sample shown in figure 2. The peak-to-peak noise in the data is consistent with the RMS noise specification at 1 s/point. This shows that increasing signal averaging improves signal-to-noise although at the expense of increased measurement time.



## NOISE AND LOW MOMENT

### Figure 4

Isothermal (IRM) and DC demagnetization (DCD) remanence curves for the 20  $\mu\text{emu}$  CoPt bit-patterned magnetic media (thin film) sample shown in figures 2 and 3. IRM is measured after the application and removal of a field with the sample initially demagnetized, and DCD is measured from the saturated state by application of increasing demagnetizing fields. Remanence curve measurement determines only the irreversible component of magnetization and thus enables the phenomena of switching to be deconvoluted from the hysteresis measurement, which generally includes a reversible component.



## NOISE AND LOW MOMENT


### Figure 5

1 h 20 min at 500 ms/point measurement of 100 first-order-reversal-curves (FORCs) for the 20  $\mu\text{emu}$  CoPt bit-patterned magnetic media (thin film) sample shown in figures 2, 3, and 4. FORCs provide information regarding magnetic interactions and coercivity distributions in magnetic materials, information that cannot be obtained from the hysteresis loop alone. In this series of FORCs there are 8,816 points. The fast (1 T/s) field sweep capability of the 8600 together with its high sensitivity makes it possible to record FORCs in a fraction of the time it takes in other magnetometer systems.

## NOISE AND LOW MOMENT

### Figure 6


10 min 40 s at 1 s/point hysteresis loop for a 20  $\mu\text{emu}$  perpendicular magnetic anisotropy (PMA) thin film. The blue curve is the measured data, which includes a diamagnetic component due to both the sample holder and sample substrate. The orange curve shows background corrected results, which removes the diamagnetic component yielding the hysteresis loop for the film only.



## NOISE AND LOW MOMENT

### Figure 7


14 min 18 s at 100 ms/point and at 5 mOe field steps for a 40  $\mu\text{emu}$  NiFe (3 nm) thin film. This demonstrates the ability of the 8600 to measure low moment magnetically soft (low coercivity) materials with high precision.




## NOISE AND LOW MOMENT

### Figure 8

28 min at 2 s/point hysteresis loop for a <2  $\mu\text{emu}$  synthetic antiferromagnetic thin film [Ta (2.5 nm)/Ru (5 nm)/Co (5 nm)/Ru (0.8 nm)/Co (5 nm)/Cu (6 nm)/Co (5 nm)/Ru (1.4 nm)/Co (10 nm)/Ta (5 nm)].






## NOISE AND LOW MOMENT

### Figure 9

Comparison of 8600 VSM (2 s/point) and AGM (100 ms/point) results for the  $<2 \mu\text{emu}$  synthetic antiferromagnetic thin film [Ta (2.5 nm)/Ru (5 nm)/Co (5 nm)/Ru (0.8 nm)/Co (5 nm)/Cu (6 nm)/Co (5 nm)/Ru (1.4 nm)/Co (10 nm)/Ta (5 nm)] shown in figure 8. Note that while the signal-to-noise are comparable in both measurements, the AGM data were recorded in 1/20th the time (1.4 min).



## NOISE AND LOW MOMENT

### Figure 10


5.5 h at 10 s/point FORC ( $N = 50$ ) for the  $<2 \mu\text{emu}$  synthetic antiferromagnetic thin film [Ta (2.5 nm)/Ru (5 nm)/Co (5 nm)/Ru (0.8 nm)/Co (5 nm)/Cu (6 nm)/Co (5 nm)/Ru (1.4 nm)/Co (10 nm)/Ta (5 nm)] shown in figures 8 and 9.

## Typical Fast Hysteresis and FORC Measurement Results

### FAST HYSTERESIS AND FORC

#### Figure 1


14 s hysteresis loop measurement (802 points, 10 ms/point) for a magnetic tape.




### FAST HYSTERESIS AND FORC

#### Figure 2

26 s hysteresis loop measurement (2000 points, 10 ms/point) for a bilayer thin film.






## FAST HYSTERESIS AND FORC

### Figure 3

3 min 30 s measurement of FORCs ( $N = 40$ , 1,376 points, 50 ms/point) for a magnetic tape. The fast (1 T/s) field sweep capability of the 8600 together with its high sensitivity makes it possible to record FORCs in a fraction of the time it takes in other magnetometer systems.



## FAST HYSTERESIS AND FORC

### Figure 4


4 min 30 s measurement of FORCs ( $N = 46$ , 2,014 points, 30 ms/point) for a magnetic stripe. The fast (1 T/s) field sweep capability of the 8600 together with its high sensitivity makes it possible to record FORCs in a fraction of the time it takes in other magnetometer systems.

## Typical Low Coercivity Measurement Results

### LOW COERCIVITY

#### Figure 1


Hysteresis loop at 1 mOe steps for a low coercivity Fe microwire. The jumps in the curve are due to Barkhausen noise, which occurs when individual magnetic domains spontaneously flip their magnetization direction.




### LOW COERCIVITY

#### Figure 2

Hysteresis loop at 5 mOe steps for a NiFe thin film.





## LOW COERCIVITY

### Figure 3


FORC (N = 100) for the NiFe thin film shown in figure 2. This shows the ability of the 8600 to measure FORCs on even very low coercivity materials.

## Typical Magnetic Multilayer Measurement Results

### MAGNETIC MULTILAYER

#### Figure 1


Hysteresis loops versus angle for a spin-valve stack. When the magnetic field is applied parallel to the exchange bias field the loop is shifted towards the right (positive field values). When applied anti-parallel the loop is shifted towards the left (negative field values). When applied perpendicular the loop is centered at  $H = 0$ .




### MAGNETIC MULTILAYER

#### Figure 2

Hysteresis loop for a [FeNi/IrMn] spin-valve. The extra steps in both the descending and ascending branches of the magnetization curve in the 3rd quadrant are related to microstructural defects/roughness of the antiferromagnetic (AFM)/ferromagnetic (FM) interfaces.






## MAGNETIC MULTILAYER

### Figure 3

FORC (N = 151) for the [FeNi/IrMn] spin-valve shown in figure 2.



## MAGNETIC MULTILAYER

### Figure 4


Hysteresis loop for a magnetic tunnel multi-junction stack.

## Typical Variable Temperature Results

### VARIABLE TEMPERATURE

#### Figure 1


Hysteresis loops at 100, 200 and 300 K (SSVT) for a core  $\text{Fe}_3\text{O}_4$  (10 nm)/shell  $\text{CoO}$  (3 nm) exchange-coupled permanent magnet. There is currently tremendous interest in developing soft (low coercivity)/hard (high coercivity) exchange-coupled magnets as they have the potential to reduce, and potentially eliminate rare-earth constituents that are used in rare-earth magnets ( $\text{NdFeB}$ ,  $\text{SmCo}$ ).




### VARIABLE TEMPERATURE

#### Figure 2

Hysteresis loops at 300, 375, 450 and 525 K (SSVT) for Fe-doped (2%)  $\text{MnCo}$  nanoparticles (NPs). Doped magnetic nanoparticles (DMNP) are a new class of ferromagnetic nanomaterials that a






## VARIABLE TEMPERATURE

### Figure 3

High temperature (oven) hysteresis loops at 25, 200, 400, 600 and 800 °C for Fe-doped (0.5%) MnCo nanoparticles (NPs). As the temperature increases the saturation magnetization ( $M_s$ ), remanence ( $M_r$ ) and coercivity ( $H_c$ ) decrease. The temperature at which a magnetic material transitions from a ferromagnetic (FM) to a paramagnetic (PM) state is called the Curie temperature  $T_{Curie}$ , this material is paramagnetic at 800 °C, which is  $>T_{Curie}$ .



## VARIABLE TEMPERATURE


### Figure 4

High temperature (oven)  $M_s(T)$  for Fe-doped (1%) MnCo nanoparticles (NPs).  $T_{Curie}$  for this material is  $\sim 925$  °C

## VARIABLE TEMPERATURE

### Figure 5


Low temperature zero-field-cooled (ZFC)/field-cooled (FC) curves for the  $\text{Fe}_3\text{O}_4$  core (10 nm)/CoO shell (3 nm) nanoparticles (NPs) shown in figure 1. The peak in the ZFC curve at each applied field corresponds to the Blocking temperature  $T_B$ . In the 8600 SSVT option with 705 gas controller and 336 temperature controller, measurements of this type are fully automated from low (80 K) to high (950 K), and high to low temperatures.



## VARIABLE TEMPERATURE

### Figure 6


FORC ( $N = 60$ ) at 120 K (SSVT) for a barium-hexaferrite nanoparticle (NP) soft/hard exchange-coupled permanent magnet. Exchange-coupled magnets can reduce, and potentially eliminate rare-earth constituents that are used in rare-earth magnets (NdFeB, SmCo).




## VARIABLE TEMPERATURE

### Figure 7

Hysteresis loop at 100 K (SSVT) for a 20  $\mu\text{emu}$  CoPt thin film. This shows that in the SSVT it is possible to measure very low moment samples with good signal-to-noise.




## Typical Rotation and Vector Measurement Results



### ROTATION AND VECTOR

**Figure 1**


$M_x(H)$  versus angle for a magnetic stripe.



### ROTATION AND VECTOR

**Figure 2**

$M_y(H)$  versus angle for the magnetic stripe shown in figure 1. The  $M_x(H)$  and  $M_y(H)$  loops at each angle are recorded simultaneously, and the angle is varied automatically (autorotation).



### ROTATION AND VECTOR


**Figure 3**

$M_y(\theta)$  versus applied field for a magnetic stripe exhibiting 2-fold anisotropy. This type of data can be used to derive the torque curve and anisotropy constant of an anisotropic material because the torque density in dyne-cm is:  $\tau(\theta) = M_y(\theta)H_x$ .

## ROTATION AND VECTOR

### Figure 4


$M_y(\theta)$  at 20 kOe for a magnetic stripe exhibiting rotational hysteresis. In a rotational hysteresis measurement the angle is swept from  $0^\circ$  to  $360^\circ$  to  $0^\circ$  under full automation (autorotation).




## ROTATION AND VECTOR

### Figure 5

$M_x(H)$  and  $M_y(H)$  for a CoPt thin film. Note that  $M_y < 10 \mu\text{emu}$  and thus demonstrates the sensitivity of the 8600 y-vector coils.






## ROTATION AND VECTOR

**Figure 6**

$M_x(H)$  versus angle for a low coercivity sputtered NiFe thin film.



## ROTATION AND VECTOR

**Figure 7**


$M_y(H)$  versus angle for the same low coercivity sputtered NiFe thin film shown in figure 6.

## MINOR LOOP

### Figure 1

Descending minor loops ( $N = 75$ ) for a magnetic stripe. Minor loops probe reversible versus irreversible magnetic processes in materials.


## Typical Minor Loop Measurement Results



## MINOR LOOP

### Figure 2

Low moment ( $30 \mu\text{emu}$ ) ascending minor loops ( $N = 12$ ) for a CoPt thin film. This data shows that minor loop measurements for low moment samples with high signal-to-noise are possible using the 8600.



## Learn more about FORC. Get our ebook.

First-order-reversal-curve (FORC) measurements provide far more information than a simple hysteresis loop. They have become indispensable for characterizing interactions and coercivity distributions in a wide array of magnetic materials.

To explore how this elegant material measurement and characterization tool can benefit your magnetism research, submit your info and receive our 36-page ebook today.

The ebook discusses FORC theory and presents measurement results for:

- Nanoparticle-based nanocomposites
- Permanent magnets
- Multi-phase ferrite magnets
- Various other materials

Go to [www.lakeshore.com/FORC](http://www.lakeshore.com/FORC) to download the ebook

### APPLICATION NOTE




## Magnetometry Measurements and First-Order Reversal Curve (FORC) Analysis

B. C. Dodrill

### First-Order-Reversal-Curve (FORC)


A FORC is measured by saturating a sample in a field  $H_{sat}$ , decreasing the field to a reversal field  $H_r$ , then sweeping the field back to  $H_{sat}$  in a series of regular field steps  $H_b$ . This process is repeated for many values of  $H_r$  yielding a series of FORCs. The measured magnetization at each step as a function of  $H_r$  and  $H_b$  gives  $M(H_r, H_b)$  which is then plotted as a function of  $H_r$  and  $H_b$  in field space. The FORC distribution  $\rho(H_r, H_b)$  is the mixed second derivative, that is,  $\rho(H_r, H_b) = -(1/2)\partial^2 M(H_r, H_b) / \partial H_r \partial H_b$ .

A FORC diagram is a 2D or 3D contour plot of  $\rho(H_r, H_b)$  with the axis rotated by changing coordinates from  $(H_r, H_b)$  to  $H_c = (H_r + H_b)/2$  and  $H_i = (H_r - H_b)/2$  where  $H_c$  corresponds to the distribution of interaction or reversal fields, and  $H_i$  the distribution of switching or coercive fields.



A typical set of FORC measurements

An example of a 3D FORC diagram



### Overview


Magnetometers are used to characterize magnetic material properties. Magnetometry techniques can be broadly classified into two categories: inductive and force-based. Common inductive methods include vibrating sample magnetometry (VSM), extraction magnetometry, AC susceptibility, and superconducting quantum interference device (SQUID) magnetometry. The two most commonly used inductive techniques are VSM and SQUID magnetometry. Alternating gradient magnetometry (AGM) is the most often used force-based technique. The measurement most commonly performed to characterize a material's magnetic properties is that of a major hysteresis loop. The hysteresis or M(H) loop is typically used to determine a material's saturation magnetization  $M_s$  (the magnetization at maximum applied field), remanence  $M_r$  (the magnetization at zero applied field after applying a saturating field), and coercivity  $H_c$  (the field required to demagnetize the material). More complex magnetization curves covering states with field and magnetization values located inside the major hysteresis loop, such as FORCs, can provide additional information that can be used to characterize magnetic interactions and coercivity distributions in magnetic materials.<sup>1</sup>

This e-book discusses some of the techniques currently used in the characterization of magnetic materials, and demonstrates the utility of the first-order-reversal-curve (FORC) measurement and analysis technique by presenting results for various magnetic materials.

## About the author

**Brad C. Dodrill** is a Senior Scientist at Lake Shore. He graduated from The Ohio State University in 1982 with a BSc degree in Physics and a minor in Mathematics. He completed 2 years of graduate studies in Physics and Electrical Engineering and took a position with Lake Shore Cryotronics in 1984 as a Research Scientist.

In his technical capacity, he is active in applications and product development initiatives in the areas of magnetic and electronic measurements and materials. He has 61 paper publications to his credit, holds 3 U.S. patents, and has lectured at numerous universities and technical conferences in the U.S., Europe, and Asia.



## About Lake Shore Cryotronics

Lake Shore Cryotronics, a DwyerOmega company, is a leading innovator in measurement and control solutions for materials characterization under extreme temperature and magnetic field conditions. High-performance product solutions from Lake Shore include cryogenic temperature sensors and instrumentation, magnetic test and measurement instruments, probe stations, and precision materials characterizations systems that explore the electronic and magnetic properties of next-generation materials. Lake Shore serves an international base of research customers at leading university, government, aerospace, and commercial research institutions and is supported by a global network of sales and service facilities.

### Contact us:

+1 614 891 2244  
sales@lakeshore.com  
www.lakeshore.com



© 2026 Lake Shore Cryotronics, LLC  
All rights reserved. Specifications are subject to change. 102726