

Instructions, Feed-Through Optic Kit Installation, DC5297

	Lake Shore Part#	Item	Qty
These components pre-assembled	DC5250	Tower, optical access	1
	P12035	O-ring, 19 mm ID x 2 mm, Buna	1
	P12036	Window, optical access	1
	DC5251	Window holder, optical	1
	P10226	Screw, M3 x 0.5 x 10, SS SHC	4
Loose parts	P10099	O-ring, -017 Buna-N	1
	P10226	Screw, M3 x 0.5 x 10, SS SHC	4

Assembly Detail: Shown upside down for clarity.

